

Product specification

Product: SALVEQUICK TEXTILE
ADHESIVE BANDAGE

Reference code: 6444

1. General description:

The dressing consists of a skin coloured elastic textile backing, evenly spread with a hypoallergenic hotmelt adhesive coating.

An absorbant pad is fixed centrally (in cross direction) to the dressing. The adhesive is covered with a release coated foil that consists of two parts that partly overlap each other. One is sealed with hotmelt on the coldsealpaper.

When taking out the dressing of the coldsealpaper, the release coated plastic foil will remain behind. One half of the dressing can now be put on the skin by one hand. The other release coated paper can be taken away. The dressing is ready for use.

The dressings are packed in a cardboard carton refill including adhesive bandages of two sizes; medium and senior.

The refill fits into the Cederroth dispenser for adhesive bandages.

2. Composition:

| | | |
|----------------|---------------------------------|-----------|
| Art # 45164400 | Medium textile adhesive bandage | 24 (encl) |
| Art # 45164401 | Senior textile adhesive bandage | 16 (encl) |

3. Technical description:

| | | |
|-----------------|---|---|
| Backing: | Elastic textile fabric Colour: Skin coloured | |
| Adhesive: | Hypoallergenic hotmelt adhesive | |
| Woundpad: | Dry laid, needle punched non woven covered with a wound contact layer | |
| Facing paper: | Singel sided siliconised release coated foil | |
| Wrapping paper: | Plain coldseal paper | |
| Dimensions: | 24 Medium | Dressing: 72 x 19 mm Wound pad: 24 x 19 mm |
| | 16 Senior | Dressing: 72 x 25 mm Wound pad: 24 x 25 mm |

4. Labelling:

| | |
|------------|---|
| On refill: | Cederroths address, telephone- and faxnumber Product name Reference code LOT-number, Expiration date and EAN-code CE mark Only use once |
| On box: | Cederroths address, telephone- and faxnumber Product name Reference code No of refills per box: 6 LOT-number and EAN-code CE mark |



On transport packaging: Cederroths address, telephone- and faxnumber
Product name
Reference code
LOT-number
EAN-code
Number of products

5. Storage conditions: The adhesive bandages should be stored in the original package for maximum protection until time of use.

6. Authorisation:

Authorized by: Yvonne Törnevik
Quality Assurance Manager

Date: 2000-09-01
