

## Product specification

---

<b>Product:</b>	CEDERROTH SMALL BLOOD-STOPPER WITH PURSE
-----------------	--

<b>Reference code:</b>	1911FR
------------------------	--------

---

<b>1. General description</b>	The idea of the small bloodstopper is that only one dressing should replace a great number of dressings, compresses and pads. In this way many hard choices can be avoided in a acute First Aid-situation.
-------------------------------	--

---

<b>2. Intended use:</b>	<p>The small bloodstopper is a universal single use First Aid-dressing with three functions:</p> <ol style="list-style-type: none"><li>1. Dressing for toes and fingers.</li><li>2. Protective dressing for superficial wounds.</li><li>3. Burns dressing with special woundprotecting surface.</li></ol>
-------------------------	---

---

<b>3. Composition:</b>	Art # 51616560      Mini bloodstopper    1
	Art # 51616570FR    Purse FR                                1

---

### 4. Labelling

On small bloodstoppers:	Cederroth, address, telephone- and faxnumber Article number of the mini bloodstopper Supplier's lotnumber CIP's lotnumber Expire date Label for 'Only use once' Sterile CE-mark
-------------------------	--

---

---

On purses:	Cederroth, address, telephone- and faxnumber Name of the product Reference code EAN-code Directions for use Intended use Composition Sterile CE-mark Article number of the purse
------------	---

On brown corrugated boxes: (package of transportation)	Cederroth, address, telephone- and faxnumber Name of the product Reference code LOT-number and EAN-code Number of products
---	--

---

## 5. Sterilization

The small bloodstoppers are sterilized by electron beam with a minimum dose of 25 kGray.

---

## 6. Storage conditions

The small bloodstoppers should be stored in the original packaging, for a maximum protection until time of use.

---

## 7. Authorisation

Authorised by:	Veronica Asp <i>Quality Assurance Manager</i>
----------------	--

Date: 981103

---