

Product specification

Product: CEDERROTH SMALL BLOOD-STOPPER WITH PURSE

Reference code: 1911FR

1. General description The idea of the small bloodstopper is that only one dressing should replace a great number of dressings, compresses and pads. In this way many hard choices can be avoided in a acute First Aid-situation.

2. Intended use: The small bloodstopper is a universal single use First Aid-dressing with three functions:

1. Dressing for toes and fingers.
 2. Protective dressing for superficial wounds.
 3. Burns dressing with special woundprotecting surface.
-

3. Composition:

Art # 51616560	Mini bloodstopper	1
Art # 51616570FR	Purse FR	1

4. Labelling

On small bloodstoppers:

- Cederroth, address, telephone- and faxnumber
- Article number of the mini bloodstopper
- Supplier's lotnumber
- CIP's lotnumber
- Expire date
- Label for 'Only use once'
- Sterile
- CE-mark



On purses: Cederroth, address, telephone- and faxnumber
Name of the product
Reference code
EAN-code
Directions for use
Intended use
Composition
Sterile
CE-mark
Article number of the purse

On brown corrugated boxes:
(package of transportation) Cederroth, address, telephone- and faxnumber
Name of the product
Reference code
LOT-number and EAN-code
Number of products

5. Sterilization

The small bloodstoppers are sterilized by electron beam with a minimum dose of 25 kGray.

6. Storage conditions

The small bloodstoppers should be stored in the original packaging, for a maximum protection until time of use.

7. Authorisation

Authorised by: Veronica Asp
Quality Assurance Manager

Date: 981103
