

## **Product specification**

Product	CEDERROTH WALL BRACKET FOR EYE WASH		
Reference code	7200		
1. General description	A holder in green hard plastics for storing a Cederroth Eye Wash bottle. The holder can be mounted on any vertical surface e.g. on a wall, in a vehicle, at the laboratory bench or on the Cederroth First Aid Panel. A Mirror is included with the Wall Bracket for visual monitoring of the rinsing.		
2. Intended Use	<ul><li>The Cederroth Wall Bracket for Eye Wash is intended for storage of a Cederroth Eye Wash bottle, to keep it organized and easy to find in an emergency situation.</li><li>The Cederroth Wall Bracket works as an opening device for the Cederroth Eye Wash.</li><li>The Cederroth Wall Bracket also works as an anti-theft system, as the Eye Wash bottle is automatically opened when removed from the holder and thus can not be returned into the holder.</li></ul>		
3. Composition	Inner package PE bagWall bracket (Polyamid plus Glassfibre)dimension (w x h x d) 84 x 262 x 70 mm (1 encl.)Mirror (styrene) with self adhesivetape on backdimension 110 x 250 x 2 mm (1 encl.)Outer transport packageNumber of inner packages pertransport package (10)		



## 4. Labelling

On Wall Bracket:	Cederroth's address, telephone and fax-number CE mark Reference code for Eye Wash refill		
On Inner package PE bag	Cederroths' address, telephone and fax-number Reference code Number of products LOT-number and EAN-code		
On Transport package:	Cederroths' address, telephone and fax-number Reference code Number of products LOT-number and EAN-code		
5. Storage conditions	Store in original packaging for maximal protection until use. No other special storing conditions.		
6. Authorisation			
	Authorized by:		
		Yvonne Törnevik Quality Assurance Manager	
	Date:	2002-07-05	