

Product specification

Product: SALVEQUICK PLASTIC
ADHESIVE BANDAGE

Reference code: 6036

1. General description:

The dressing consists of a skin coloured PE (polyethylene) backing, evenly spread with a polyacrylate adhesive coating.

An absorbant pad is fixed centrally to the dressing.

The adhesive is covered with a release coated foil that consists of two parts that partly overlap each other. One is sealed with hotmelt on the coldsealpaper.

When taking out the dressing of the coldsealpaper, the release coated plastic foil will remain behind. One half of the dressing can now be put on the skin by one hand. Then the other release coated paper can be taken away. The dressing is ready for use.

The dressings are packed in a cardboard carton refill including adhesive bandages of two sizes; medium and senior.

The refill fits into the Cederroth dispenser for adhesive bandages.

2. Composition:

Art # 45160330	Medium plastic adhesive bandage	27 (encl)
Art # 45160331	Senior plastic adhesive bandage	18 (encl)

3. Technical description:

Backing:	Polyethylene foil. Colour: Skin coloured	
Adhesive:	Transparent polyacrylate adhesive	
Woundpad:	Dry laid, needle punched non woven covered with a wound contact layer	
Facing paper:	Singel sided siliconised release coated paper/foil	
Wrapping paper:	Plain coldseal paper	
Dimensions:	27 Medium	Dressing: 72 x 19 mm Wound pad: 24 x 19 mm
	18 Senior	Dressing: 72 x 25 mm Wound pad: 24 x 25 mm

4. Labelling:

On refill:	Cederroths address, telephone- and faxnumber Product name Reference code LOT-number, Expiration date and EAN-code CE mark Only use once
On box:	Cederroths address, telephone- and faxnumber Product name Reference code No of refills per box: 6 LOT-number and EAN-code CE mark



On transport packaging: Cederroths address, telephone- and faxnumber
Product name
Reference code
LOT-number
EAN-code
Number of products

5. Storage conditions: The adhesive bandages should be stored in the original package for maximum protection until time of use.

6. Authorisation:

Authorized by: Yvonne Törnevik
Quality Assurance Manager

Date: 2004-03-24
